

LIQUIPETTE®

Plastic Transfer Pipettes

Manufactured from virgin LDPE, free from phthalate and silicone

Our Liquipette® plastic disposable Pasteur pipettes have a well-earned reputation for quality and reliability.



Non-sterile versions P511, P503, P406 and P509 are supplied in cartons. The unique press-out section enables easy access and keeps the contents dust and particle free when not in use.

Available in 14 styles and 25 packaging options including fine points, long-lengths, small bulbs, and special styles for use with diagnostic kits.

- Economical and safe alternative to glass Pasteur pipettes
- Suitable for contact with foodstuffs
- Manufactured in facilities operating to ISO 9001-2015 and ISO 13485 standards
- Sterile and non-sterile versions, in a variety of packaging options



Liquipette® pipettes are available worldwide with special pricing available for bulk orders; please enquire for more information.

LIQUIPETTE® Plastic Transfer Pipettes


Product	P522 LIQUIPETTE® Mini, Graduated to 0.5ml	P423 LIQUIPETTE® Fine Point, Micro Extended	P523 LIQUIPETTE® Mini-Plus, Graduated to 0.5ml	P520 LIQUIPETTE® Junior, with Mixer	P521 LIQUIPETTE® Junior, Fine Point, with Mixer	P409 LIQUIPETTE® Fine Point, Small Bulb	P509 LIQUIPETTE® 1ml Graduated, Small Bulb	P406 LIQUIPETTE® Fine Point	P507 LIQUIPETTE® Thin Stem		
Length (cm)	8.5	10.4	11.2	12.2	13.0	13.2	14.0	15.0	15.2		
Bulb Draw (ml)	1.0	1.0	1.0	0.5	0.5	1.5	3.0	3.5	3.5		
Drop Size* (µl)	40	20	40	50	25	20	40	20	40		
PACKAGING	NON-STERILE	Cat No Unit Qty	127-P522-000 1000	127-P423-000 400	127-P523-000 1000	127-P520-000 1000	127-P521-000 1000	127-P409-000 500	127-P509-000 500	127-P406-000 500	127-P507-000 500
	STERILE										
	SINGLE WRAP	Cat No Unit Qty		127-P423-STR 250				127-P509-STR 250			
	10 PACK	Cat No Unit Qty							127-P406-10S 1000		
	20 PACK	Cat No Unit Qty		127-P423-20S 500					127-P406-20S 1000		

* Drop sizes are nominal and may vary according to operator technique.
Sterile products are γ irradiated.

= ACTUAL SIZE

Liquipette® is a registered trademark of Elkay Laboratory Products (UK) Ltd.
Drawings are shown for illustration purposes only.
We reserve the right to alter product dimensions accordingly.

Product	P511 LIQUIPETTE® 1ml Graduated	P503 LIQUIPETTE® 3ml Graduated	P190 LIQUIPETTE® 10ml Large Bulb, Thin Stem	P514 LIQUIPETTE® Extra-Long Stem	P330 LIQUIPETTE® Maxi-30, 10ml Draw
Length (cm)	15.2	15.2	17.0	22.9	30.0
Bulb Draw (ml)	3.5	3.5	10.0	3.5	10.0
Drop Size* (µl)	45	45	55	40	40



PACKAGING	NON-STERILE	Cat No Unit Qty	127-P511-000 500	127-P503-000 500	127-P190-000 250	127-P514-000 500	127-P330-000 100	
	STERILE							
	SINGLE WRAP	Cat No Unit Qty	127-P511-STR 250	127-P503-STR 250		127-P514-STR 250	127-P330-STR 150	
	10 PACK	Cat No Unit Qty	127-P511-10S 1000	127-P503-10S 1000				
	20 PACK	Cat No Unit Qty	127-P511-20S 1000	127-P503-20S 1000				

Elkay c/o Calibre Scientific,

Unit 11 A&B, Telford Way (off Whittle Way), The Advanced Manufacturing Park, Catcliffe, Rotherham, S60 5GP